



















**EN**  
**Symbols Glossary**

The following symbols may be used in product labeling of Oertli Instrumente AG:

Symbol	Description of symbol	Standard Reference
	<b>Manufacturer</b> Indicates the medical device manufacturer	ISO 15223-1:2021 Reference number of symbol: 5.1.1
	<b>Authorized representative in the European Community/ European Union</b> Indicates the authorized representative in the European Community/ European Union	ISO 15223-1:2021 Reference number of symbol: 5.1.2
	<b>Use-by date</b> Indicates the date after which the medical device is not to be used.	ISO 15223-1:2021 Reference number of symbol: 5.1.4
	<b>Batch code</b> Indicates the manufacturer's batch code so that the batch or lot can be identified.	ISO 15223-1:2021 Reference number of symbol: 5.1.5
	<b>Catalogue number</b> Indicates the manufacturer's catalogue number so that the medical device can be identified.	ISO 15223-1:2021 Reference number of symbol: 5.1.6
	<b>Serial number</b> Indicates the manufacturer's serial number so that a specific medical device can be identified.	ISO 15223-1:2021 Reference number of symbol: 5.1.7
	<b>Sterilized using ethylene oxide</b> Indicates a medical device that has been sterilized using ethylene oxide.	ISO 15223-1:2021 Reference number of symbol: 5.2.3
	<b>Do not resterilize</b> Indicates a medical device that is not to be resterilized.	ISO 15223-1:2021 Reference number of symbol: 5.2.6
	<b>Do not use if package is damaged</b> Indicates that a medical device that should not be used if the package has been damaged or opened and that the user should consult the instructions for use for additional information.	ISO 15223-1:2021 Reference number of symbol: 5.2.8

	<p><b>Single sterile barrier system</b> Indicates a single sterile barrier system</p>	<p>ISO 15223-1:2021 Reference number of symbol: 5.2.11</p>
	<p><b>Temperature Limit</b> Indicates the temperature limits to which the medical device can be safely exposed.</p>	<p>ISO 15223-1:2021 Reference number of symbol: 5.3.7</p>
	<p><b>Do not re-use</b> Indicates a medical device that is intended for one single use only.</p>	<p>ISO 15223-1:2021 Reference number of symbol: 5.4.2</p>
	<p><b>Consult instructions for use</b> Indicates the need for the user to consult the instructions for use.</p>	<p>ISO 15223-1:2021 Reference number of symbol: 5.4.3</p>
	<p><b>Caution</b> Indicates that caution is necessary when operating the device or control close to where the symbol is placed, or that the current situation needs operator awareness or operator action in order to avoid undesirable consequences.</p>	<p>ISO 15223-1:2021 Reference number of symbol: 5.4.4</p>
	<p><b>Medical device</b> Indicates the item is a medical device.</p>	<p>ISO 15223-1:2021 Reference number of symbol: 5.7.7</p>
	<p><b>Unique device identifier</b> Indicates a carrier that contains unique device identifier information.</p>	<p>ISO 15223-1:2021 Reference number of symbol: 5.7.10</p>
	<p><b>Keep away from sunlight</b> Indicates a medical device that needs protection from light sources</p>	<p>ISO 15223-1:2021 Reference number of symbol: 5.3.2</p>
	<p><b>Observe instructions for use</b></p>	<p>IEC 60601-1:2005/A1:2012, Table D.2, No. 10</p>



**Oertli Instrumente AG**  
Hafnerwisenstrasse 4  
9442 Berneck  
Schweiz

