












## EN Symbols Glossary

The following symbols may be used in product labeling of Oertli Instrumente AG:

Symbol	Description of symbol	Standard Reference
	<b>Manufacturer</b> Indicates the medical device manufacturer, as defined in EU Directive 93/42/EEC	ISO 15223-1:2016 Reference number of symbol: 5.1.1
	<b>Use-by date</b> Indicates the date after which the medical device is not to be used.	ISO 15223-1:2016 Reference number of symbol: 5.1.4
	<b>Batch code</b> Indicates the manufacturer's lot number so that the batch or lot can be identified.	ISO 15223-1:2016 Reference number of symbol: 5.1.5
	<b>Catalogue number</b> Indicates the manufacturer's catalogue number so that the medical device can be identified.	ISO 15223-1:2016 Reference number of symbol: 5.1.6
	<b>Serial number</b> Indicates the manufacturer's serial number so that a specific medical device can be identified.	ISO 15223-1:2016 Reference number of symbol: 5.1.7
	<b>Sterilized using ethylene oxide</b> Indicates a medical device that has been sterilized using ethylene oxide.	ISO 15223-1:2016 Reference number of symbol: 5.2.3
	<b>Do not use if package is damaged</b> Indicates a medical device that should not be used if the package has been damaged or opened.	ISO 15223-1:2016 Reference number of symbol: 5.2.8
	<b>Do not re-use</b> Indicates a medical device that is intended for one use, or for use on a single patient during a single procedure.	ISO 15223-1:2016 Reference number of symbol: 5.4.2
	<b>Consult instructions for use</b> Indicates the need for the user to consult the instructions for use.	ISO 15223-1:2016 Reference number of symbol: 5.4.3
	<b>Caution</b> Indicates the need for the user to consult the instructions for use for important cautionary information such as warnings and precautions that cannot, for a variety of reasons, be presented on the medical device itself.	ISO 15223-1:2016 Reference number of symbol: 5.4.4

	<b>Observe instructions for use</b>	IEC 60601-1:2005/A1:2012, Table D.2, No. 10
---	-------------------------------------	--



**Oertli Instrumente AG**  
Hafnerwisenstrasse 4  
CH-9442 Berneck  
Schweiz