



Endo Illuminators

MAKING THE DIFFERENCE WITH ENHANCED INTRA- OCULAR ILLUMINATION




Endo Illuminators



Features and benefits

- High light output
- Enhanced stiffness of 25 G light instruments for better manipulation
- Panoramic tip for wide field illumination
- Shielded design to avoid glare
- Comfortable and sterile connection with Oertli comfort-connector
- Color coded by gauge size

Two types of illuminators for enhanced intra-ocular illumination



**Endo Illuminator 90°,
comfort-connector**
VV300101


⊗ ● 20 G

**Endo Illuminator 90°,
comfort-connector**
VV300103

⊗ ● 23 G

**Endo Illuminator 90°,
comfort-connector**
VV300105

⊗ ● 25 G



**Endo Illuminator Panorama, shielded,
comfort-connector**
VV300191

⊗ ● 20 G

**Endo Illuminator Panorama, shielded,
comfort-connector**
VV300194

⊗ ● 23 G

**Endo Illuminator Panorama, shielded,
comfort-connector**
VV300195

⊗ ● 25 G

All articles are delivered sterile, single use, box of 10 and are also available in various vitrectomy packs

Stiff and more rigid 25 G endo illuminators

The 25 G light instruments feature a stiffer shaft design. An additional cannula increases the rigidity of the light pipe and prevents the fiber from bending when manipulating.




New comfort-connector for easy handling

The new comfort-connector maximizes light output and allows an easy installation out of the sterile field.



Safe protection and easy preparation


All illuminators will come in a protective packaging. This will facilitate the preparation for the nurse.



Compatibility


The new endo illuminators fit to all existing Oertli light output of the GoodLight LED light sources: OS4 and Faros anterior / posterior


Chandelier Illumination



Features and Benefits

- Extrawide and homogeneous illumination
- Accurate positioning with bendable loop of the fiber
- Reliable snap-lock into all Oertli Caliburn trocars
- Compatible to Oertli Caliburn 23 G and 25 G trocars
- Precise and comfortable connector to Oertli surgery platforms
- Postoperative tight incision



Chandelier light with trocar,  ● 25 G
comfort connector
VV300152

All articles are delivered sterile, single use, box of 10

oertli[®]
S W I T Z E R L A N D

Oertli Instrumente AG

Hafnerwisenstrasse 4
9442 Berneck
Switzerland

T +41 71 747 42 00
F +41 71 747 42 90

www.oertli-instruments.com

Oertli[®] and the Oertli logo are registered trademarks of Oertli Instrumente AG.
Caliburn[™] and GoodLight LED[™] are trademarks of Oertli Instrumente AG.

Not available for sales in the US