


Intraocular illumination

MAKING THE DIFFERENCE WITH TROCAR BASED CHANDELIER VITRECTOMY



NEW

One-step access to
bimanual surgery


Flexible chandelier


Using chandelier illumination provides a true bimanual manipulation during vitreoretinal surgery. Our trocar based solution offers a one-step access for easy and quick connection.


There are many different surgical situations, in which chandelier illumination is beneficial for improving intraocular visibility and thereby achieving favourable surgical outcomes. It facilitates clear intraocular visualization and the surgeon can perform manipulations bimanually and efficiently.

Features and benefits for all Oertli chandeliers


- One-step access with quick trocar based solution
- Accurate fixation due to reliable snap lock connection into all Oertli Caliburn trocars
- Extrawide and homogeneous illumination
- Postoperative tight sclerotomy
- Precise and comfortable magnetic connection to Oertli surgery platforms



Chandelier illuminator flexible,  ● 25 G
with Caliburn trocar, snap lock
for 23G and 25G
VV300155

Chandelier illuminator flexible,  ● 27 G
with Caliburn trocar, snap lock
for 23 G, 25 G, 27 G
VV300157

All articles are delivered sterile, single use, box of 10





Features and Benefits



- Global visibility of the fundus for all surgical maneuvers
- High flexibility in illumination through flexible loop of fiber
- Compatible to all Oertli 23G, 25G and 27G trocars

Scan the QR code
to find out more.



Bendable chandelier



Chandelier illuminator bendable,   25 G
with Calibur trocar, snap lock
for 23 G and 25 G
VV300152

All articles are delivered sterile, single use, box of 10


Features and Benefits


- Rigid positioning through bendable loop of fiber
- Flexible re-positioning in all Oertli trocars
- Compatible to Oertli 23 G and 25 G trocars


Endo illuminators


Features and benefits


- Enhanced stiffness of 25 G light instruments for better manipulation
- Panoramic tip for wide field illumination
- Shielded design to avoid glare
- Comfortable and sterile connection with Oertli comfort-connector
- Each illuminator is accompanied by the unique ViPer illuminated scleral indenter





Endo illuminator 90°,  ● 20 G
with ViPer scleral indenter,
comfort-connector
VV300101


Endo illuminator panorama,  ● 20 G
shielded, with ViPer scleral
indenter, comfort-connector
VV300191

Endo illuminator 90°,  ● 23 G
with ViPer scleral indenter,
comfort-connector
VV300103

Endo illuminator panorama,  ● 23 G
shielded, with ViPer scleral
indenter, comfort-connector
VV300194

Endo illuminator 90°,  ● 25 G
with ViPer scleral indenter,
comfort-connector
VV300105

Endo illuminator panorama,  ● 25 G
shielded, with ViPer scleral
indenter, comfort-connector
VV300195

Endo illuminator widefield,  ● 27 G
shielded, with ViPer scleral
indenter, comfort-connector
VV300197

All articles are delivered sterile, single use, box of 10
and are also available in various vitrectomy packs

oertli[®]
S W I T Z E R L A N D

Oertli Instrumente AG

Hafnerwisenstrasse 4
9442 Berneck
Switzerland

T +41 71 747 42 00
F +41 71 747 42 90

www.oertli-instruments.com

Oertli® and the Oertli logo are registered trademarks of Oertli Instrumente AG.
ViPer®, Calibur™ and GoodLight LED™ are trademarks of Oertli Instrumente AG.

Not available for sales in the US